

TECHNICKY LIST / TECHNICAL SHEET / TECHNISCHES DATENBLATT

Drafted : Ing. Pavel Lubina
Approved : Jaroslav Kučienksy



Document :
CliniRoll/23/5812

Revision : rev00/23

I.

Product :

Clini – Roll is white paper in a roll, material 100% cellulose.

II.

Usage:

Clini – Roll is used on the examination bed to hygienically protect the patient during examination by the doctor.

Instructions: the medical facility worker puts the roll of paper into the holder by the examination bed, and stretches the paper over the entire area of the bed. After the examination, the paper is disposed of in sorted waste, according to the waste management of each workplace.

III.

Material:

100% cellulose, 2 ply, white

34 gr/m² layer,

Width: 39/50/59 cm and length 50 m, 20 cm length of one piece.

Each roll is hygienically wrapped in PE foil and provided with a label.

Roll diameter: 11.5 cm

IV.

Print the date on the label:

Name, REF number, width and length of the product.

Production number, Production lot.

The product is NON-STERILE,



V.

Storage

Clini Roll must be stored in a dry storage area, at a temperature of +5 / +35 C, properly ventilated and absolutely protected against moisture, rodents, and fire. Storage must be properly secured in accordance with fire and occupational safety regulations.

VI.

Reference numbers

REF	Product
CLX39	Clini Roll width 39 cm, length 50 m
CLX50	Clini Roll width 50 cm , length 50 m
CLX59	Clini Roll width 59 cm , length 50 m



VII.

Logistics Data

	39 cm	50 cm	59 cm
Colour	white	white	white
Embossing	Yes	Yes	Yes
Width 1 pcs	20 cm	20 cm	20 cm
Width roll	39 cm	50 cm	59 cm
Length 1 roll	50 m	50 m	50 m
Weight gr/m ²	34	34	34
Rolls in carton	6 ks	9 ks	9 ks
Weight 1 role	800 gr	850 gr	850 gr
Weight 1 carton	6,5 kg	7,8 kg	8,8 kg
cartons on pallet	50	20	20
Váha 1 palety	190 kg	155 kg	170 kg



IX.

Conclusion:

Clini Roll paper in a roll for examination beds meets all parameters for use in medical facilities and for general hygiene and cleaning purposes.

X.

Contact the manufacturer and responsible person:

DispoPharma s.r.o.

Rumunská 12, 120 00 Prague 2

Czech republic

ID : 05009049

VAT ID : CZ05009049

EORI customs registration : CZ05009049

E-mail : info@dispopharm.cz

phone : +420 602 243 441

Responsible person : p. Jiří Němeček



XI.



Sustainability

The DispoPharm company pays close attention to the company's sustainability and environmental protection development plan, where packaging materials are used from recyclable material and conducts regular employee training and client presentation training on proper waste management as part of sustainable development.

XII.

Sustainable Development - vision





XIII.

General use :





Headquarters :

DispoPharm s.r.o.

Rumunská 12, čp. 1720

120 00 Prague 2

Czech republic

IČ : 050 09 049

Bank connecting :

ČSOB a.s.

EUR account : 298319445/0300

CZK account : 298294485/0300

IBAN : CZ03 0300 0000 0002 9831 9945

Swift: CEKOCZPP

e-mail : info@dispopharm.cz

phone : +420 602 243 441

